



EN-E-09-003-001

Automatic disposable plastic tip loading machine



Automatize the complex operations



Small and multi-functional device

ATS-3000

New Release



 **MedicaTec Corporation**

MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,
Beijing, China, 101100

TEL: +86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicattec.cn>

E-mail: zhangxiang@medicattec.cn

ATS-3000

Automatize the complex operations

Small desk-top automatic plastic tip loading machine

Product performance

• Treating rate

Ten boxes (can be added halfway) can be placed at once, and 90 seconds will be needed for loading a whole box.

• Treatment of plastic tip

Remove the superposed tips and send to the loading site by multiple detections.

A total of 3000 200 μ l (tips of other volumes can also be produced in forms of OEM) plastic tips can be loaded at once.

Load the plastic tips from the open site,

and protect the surface membrane on plastic needle.

• Control program

Lines or rows needed to be loaded can be reset,

and 9 user procedures can be stored

- 1.Transfer box for plastic tips
- 2.Receipt box of plastic tips
- 3.Outlet of plastic tips
- 4.Transfer platform of plastic tips shelf



Product specifications

TYPE	ATS-3000
External size	655 (W) × 525 (D) × 715mm (H)
Weight	65kg
Power	AC220V, 80VA

MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901,
Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicattec.cn>

E-mail: zhangxiang@medicatec.cn