




AMW-96SX ELISA 96ch Micro plate washer with plate hopper feeder

◆◆◆◆◆◆◆◆◆◆ 96ch Metal Nozzle Liquid head ◆◆◆◆◆◆◆◆◆◆

AMW-96SX

New Release



MedicaTec Corporation

MedicaTec Corporation

Tong Zhou-Qu, Ba Li Qiao-Nan Jie No. 68, E#2901,
Beijing, China 101100

TEL:+86-10-89506062 Cell:+86-15110235386

URL:<http://www.medicattec.cn>


E-mail:zhangxiang@medicattec.cn

AMW-96SX

96ch metal double nozzle, Table top ELISA washer

【Characteristic】

- 96ch washing at same time
- Washing liquid head and drain piping automatic cleaning
- Over flow washing
- Can set the intensity of washing dispense
- Dispense, suction sensor (optional)
- Speed 20sec./plate (depend on user program)
- (Maximum 20)plates hopper (can add plates on any time on the process)



96ch metal nozzle liquid head

Optional Unit

Dispense and suction liquid sensor

【Specifications】

TYPE	AMW-96SX
Ambient	10 ~ +50℃ (None Dewfall) , humidity less than 85%
Dispense volume	50μL~990μL
Wash time	1-9time
Wash method	96ch dispense at same time, suction at same time
Microplate	All sorts of 96ch microplate
Residual liquid	Every well≤1.5μL
Washing liquid dunk time	maximum 999 second
Washing mode	Normal, Dispense, Aspirate, Dispense End, 30 user program memory
Micro plate hopper feeder	20plates supply and 20plates collect hopper (can add plates on any time on the process)
Dimension	440 (W) ×490 (D) ×260mm (H)
Weight	27Kg
Power	AC100V(AC220 available), 250VA
Standard	60L/min Vacuum pump, 5L drain bottle, 2L trap bottle

MedicaTec Corporation

Tong Zhou-Qu Ba Li Qiao-Nan Jie No. 68, E#2901, Beijing, China, 101100

TEL: + 86-10-89506062 Cell: +86-15110235386

URL: <http://www.medicathec.cn>

E-mail: zhangxiang@medicatec.cn